

STANDARD OPERATING PROCEDURE

Title: Preparation of Soluble Lysates from Tumor Tissue

SOP#: WU-SOP-TL1-01

Version #: 2

**Authors: Petra Erdmann-Gilmore
Sherri Davies, PhD**

Reviewer: Matthew Meyer, PhD

**WU-PCC Directors: Matthew J. Ellis, MB, PhD
R. Reid Townsend, MD, PhD**

Date Approved: 6/22/2014 Date Modified: 11/13/10/14

1. PURPOSE

The purpose of this document is to describe the preparation of a soluble lysate from cryopulverized breast cancer xenograft tissue in a non-denaturing lysis buffer.

2. SCOPE

This procedure is used to prepare a pooled tumor lysate from cryopulverized tissue that is used as a test matrix for the development of PRM assays for peptides quantification in tumor digests.

3. RESPONSIBILITIES

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

4. EQUIPMENT

Covaris S220X Focused Acoustics System with cooling unit
Micro-centrifuge (Eppendorf, Model No. 5424)
Rainin™ Pipet-lite XLS, P20, P200, P1000
Axygen® MAXYmum™ recovery tips;
P200 and P20: T-200-C-L-STK,
P1000: T-1000-C-L-R

5. MATERIALS

Item	Vendor	Catalog #	Unit
12 x12 glass tubes and caps	Covaris	520080	100 per pack
1.7mL Eppendorf tubes	Axygen	MCCT-175-C	
EDTA	Sigma	E1644	100g
EGTA	Sigma	E4378	25g
1M HEPES	Sigma	3537	100mL
Hydrochloric acid	Fisher	A144S-500	500 mL
Phenyl methylsulfonyl fluoride	Thermo Scientific	36978	5g
Protease Inhibitor Cocktail-Complete	Roche	11697498001	20 Tablets
Phosphatase Inhibitor Cocktail 2	Sigma	P5726	
Phosphatase Inhibitor Cocktail 3	Sigma	P0044	
Sodium Chloride	Fisher	S271	500g
Sodium Orthovanadate	Sigma	S6508	50g
Sodium Fluoride	Sigma	S-7920	100g
Triton X-100	Sigma	93443	100mL
18.2 mΩ water	Millipore		
Dry ice			
Ice			
Pipette tips	Axygen	T-1000-C-L-R	

6. REAGENTS

- A.** 100mM PMSF in ethanol.
- 1) Weigh 17.4mg of PMSF into a 1.7mL Eppendorf tube.
 - 2) Solubilize in 1 mL of 100% ethanol.
 - 3) Aliquot and store at -20°C.
- B.** HEPES, pH 7.5 (1M)
- C.** NaCl (5M)
- D.** Triton X-100 (10%)
- E.** EDTA (100 mM)
- F.** EGTA (250 mM)
- G.** NaF (100 mM)
- H.** NaVO₄ (100 mM)
- I.** Protease Inhibitor Cocktail (Roche)
- J.** Preparation of Phosphatase Inhibitor Cocktail 3
- 1) Add 500 μL to 50 mL
 - 2) Mix well and aliquot (volume) for storage at -20°C
- K.** Preparation of MIB Lysis Buffer (Prepare fresh prior to use)
- 1) 50mM HEPES, pH 7.5
 - 2) 150 mM NaCl
 - 3) 0.5% Triton X-100
 - 4) 1mM EDTA
 - 5) 1mM EGTA
 - 6) 10mM NaF

- 7) Protease inhibitor cocktail (2 tablets/50 mL)
- 8) Phosphatase inhibitor cocktail 2 (500 μ L/50 mL)
- 9) Phosphatase inhibitor cocktail 3 (500 μ L/50 mL)

Albumin Standard (Pierce. 23209)

Advanced Protein Assay Reagent (Cytoskeleton Inc. ADV01)

7. PROCEDURE

A. Set up Covaris S220X

- 1) Set up Covaris S220X. Add approximately 1.55L of DI water to the reservoir. Place reservoir with water under the transducer assembly and push transducer assembly down into the reservoir.
- 2) Check water level with 12 x 12 glass tube; the entire glass tube up to the cap should be immersed in water.
- 3) Switch on the S220X unit.
- 4) Start Sonolab 7.0 software on the laptop. The degas pump will start automatically.
- 5) Start the water chiller and set temperature to 4^oC.
- 6) The Sonolab software will show green check marks when the instrument is ready in the reservoir is complete and that the operating temperature has been reached. It will take approximately 30 minutes.

B. Protein Extraction From Tissue Sample

- 1) Collect tissue samples from -80^oC freezer and put them on dry ice until buffer is added. For the preparation of the P7 matrix, equal portions of tissue from different subtypes of breast cancer tumors were pooled into a homogenous mix prior to solubilization.
- 2) Add 400 μ L of prepared MIB lysis buffer to frozen tissue in cryovial. Transfer buffer and tissue to 12x12 glass tube with pipette. It helps to cut off a portion of the pipette tip to increase the diameter of the pipet tip to allow the tissue to be aspirated with the solubilization buffer.
- 3) Add an additional 400 μ L MIB lysis buffer to the cryovial and rinse the vial. Transfer buffer and any remaining tissue to the 12x12 glass tube.
- 4) Cap the 12x12 glass tube and place into the tube holder for the Covaris S220X.
- 5) Start extraction method (see Table 1 for method parameters).
- 6) Once extraction run is completed, remove the 12 x 12 glass tube from the S220X and transfer sample into a 1.7mL Eppendorf tube.
- 7) Spin sample at 16,000 rcf for 10 minutes.
- 8) Remove supernatant and add to a new 1.7 mL Eppendorf tube.
- 9) Sample is ready for further processing or can be stored at -80^oC.

C. Shut down Covaris S220X

- 1) When all samples are processed, shut down degas pump in the Sonolab software

- 2) Lift transducer assembly out of the water bath and empty water from reservoir.
- 3) Start the degas pump. The pump will pump out any remaining water in the degas coils.
- 4) Once the degas pump stops, start it one more time.
- 5) Lift transducer assembly out of the reservoir and remove all water. Dab transducer assembly dry with a kimwipe. Dry the reservoir with a kimwipe.
- 6) Place reservoir back under transducer assembly and lower transducer assembly.
- 7) Exit from Sonolab software, and then shut off the S220X unit and the water chiller.

D. Protein Determination

- 1) Prepare serial dilutions of 2mg/mL BSA standard in the range of 0.007 mg/mL to 2mg/mL BSA in DI water.
- 2) Dilute samples 1:10 in water; make 50 μ L.
- 3) Dilute lysis buffer 1:10 in water; make 1mL
- 4) Prepare a 1:50 dilution of each sample using the 1:10 diluted samples by pipetting 24 μ L of 1:10 diluted lysis buffer in a 0.5mL PCR tube and adding 6 μ L of the 1:10 diluted sample.
- 5) Prepare a 1:100 dilution of each sample using the 1:10 diluted samples by pipetting 27 μ L of 1:10 diluted lysis buffer in a 0.5mL PCR tube and adding 3 μ L of the 1:10 diluted sample.
- 6) Pipet 10 μ L of 1:10 diluted lysis buffer into each standard cuvette, then add 10 μ L of standard. To the blank cuvette, add 10 μ L of DI water.
- 7) Pipet 10 μ L of water into each sample cuvette, then add 10 μ L of each sample dilution. Samples are assayed in duplicate.
- 8) Add 1mL of Advanced Protein Assay Reagent, mix and read in Thermo BioMate 3 spectrophotometer at 590nm.
- 9) Plot the response curve in Excel and use the linear regression line to calculate sample concentration.

8. REFERENCED DOCUMENTS

Mertins P (2014) "Ischemia in tumors induces early and sustained phosphorylation changes in kinase pathways but does not affect global protein levels. Mol. Cellular Proteomics. 13, 1690-1704.