

# STANDARD OPERATING PROCEDURE

**Title: MRM mass spectrometry for the analysis of immuno-MRM assay samples using a 6500 QTRAP**

**Version #: 1**

**Author: Paulovich Lab**

**Date: 7/1/2016**

## Purpose

The purpose of this document is to describe a multiple reaction monitoring (MRM) method for quantitative analysis of peptides from immunoaffinity enriched samples (immuno-MRM samples). The method employs retention time scheduling for enhanced quantification.

## Scope

This procedure includes the setup of MRM methods on a 6500 QTRAP. Liquid chromatography (LC) parameters and methods are described in a separate SOP.

## Responsibilities

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

## Equipment

- Mass spectrometer: Sciex 6500 QTRAP
- Ion Source: Sciex NanoSpray III
- Emitter tip: TaperTip 20  $\mu$ m ID (NewObjective)

## Materials

- Water, Optima® LC/MS, suitable for UHPLC-UV (W6-4, Fisher Scientific)
- Acetic acid (ACS reagent,  $\geq 99.7\%$ , 242853, Sigma)
- Acetonitrile, Optima® LC/MS, suitable for UHPLC-UV (A955-4, Fisher Scientific)
- Formic acid (FA) (1.11670.1000, EMD Millipore)

## Procedure

1. Setup MRM method parameters
  - a. Source/Gas Parameters:
    - i. Curtain Gas (CUR): 20
    - ii. IonSpray Voltage (IS): 3.5 kV
    - iii. Ion Source Gas 1 (GS1): 18
    - iv. Ion Source Gas 2 (GS2): 0
    - v. Interface Heater Temperature (IHT): 160 °C
    - vi. Collision Gas (CAD): Medium
  - b. Scheduled MRM Parameters
 

Note: These parameters were used with the aim of obtaining at least 10 data points per peak.

    - i. MRM detection window: 60 sec
    - ii. Target Scan Time: 0.5 sec
  - c. MS Parameters:
    - i. Declustering Potential (DP): 100 V
    - ii. Entrance Potential (EP): 10 V
    - iii. Collision Energy (CE): From Skyline<sup>1</sup>, based on optimized values from synthetic peptides.
    - iv. Collision Cell Exit Potential Q1 (CXP): 10 V
  - d. Advanced MS Parameters:
    - i. Resolution Q1: Unit
    - ii. Resolution Q3: Unit
    - iii. Intensity Threshold (total count): 0
    - iv. Settling time: 0 ms
    - v. Pause between mass ranges: 5.007 ms
2. Test the system suitability with an appropriate standard once the column is conditioned.
 

Note: Pierce iRT standard (88320, Pierce) was injected on a conditioned column at least twice to be able to assess retention time reproducibility and whether peak shapes and intensities are acceptable based on historic data.

  - a. Retention time shift of <0.5 minutes (a third injection might be needed if the first two runs' retention times shift by >0.5 min; the LC system needs to be troubleshot if the retention times still shift >0.5 min after a third run).
  - b. Minimal tailing or fronting (check the LC connections if necessary).
  - c. No drop-out of electrospray (sparge the LC solvents or replace the emitter tip if needed).
  - d. If the peak intensities are unacceptable, troubleshoot the autosampler or clean the MS.
3. Identify scheduling times for target peptides/transitions
  - a. Target LC-MRM method preparation
    - i. Load the Skyline file containing peptides and transitions that will be monitored during the analysis.

- ii. Export a scheduled transition list from Skyline using the “values from a single data set” and choosing the correct data file.
- iii. Import the scheduled transition list into an MRM acquisition method on the 6500 with all parameters set as above (step 1)
- iv. Set up the autosampler and LC methods as in the LC SOP.
- v. Inject one immuno-MRM sample.
- b. Instrument performance evaluation and adjustment of scheduled retention times
  - i. Import the latest data file into the Skyline file.
  - ii. Check the automatic integration of all peaks.
    1. Make sure integration start and stop is identical for all transitions of a precursor (go to “Settings”, and check “Integrate All” to enable this feature automatically).
    2. Manually adjust the integration of peaks if necessary.
  - iii. Check that peak shapes are acceptable.
    1. Minimal tailing or fronting.
    2. No drop-out of electrospray (for this, no smoothing should be applied in Skyline (View/Transform/None)).
    3. No missing transitions.
  - iv. Export a new scheduled transition list from Skyline if needed (the retention times can shift based on whether peptides are tested in a buffer background or whether the matrix is more complex, such as in the eluted samples after an immunoaffinity capture from a plasma matrix).
4. Once the instrument is performing acceptably, run the remaining samples of the immuno-MRM experiments in the order given in the SOPs for response curves, validation samples, or unknown samples. The file naming convention follows:
  - a. For the response curve experiments:  
AssayGroup\_TypeOfExperiment\_AnalyteAmountLevel\_CaptureReplicate#\_InjectionReplicate#
  - b. For the repeatability (validation) experiments:  
AssayGroup\_AnalyteLevel\_Day#\_CaptureReplicate#\_InjectionReplicate#

## References

1. MacLean, B., Tomazela, D. M., Shulman, N., Chambers, M., Finney, G. L., Frewen, B., Kern, R., Tabb, D. L., Liebler, D. C., and MacCoss, M. J. (2010) Skyline: an open source document editor for creating and analyzing targeted proteomics experiments. *Bioinformatics* 26, 966-968.