

# STANDARD OPERATING PROCEDURE

**Title: Liquid Chromatography, Eksigent Ultra 2D+, trap and elute**

**Version #: 1**

**Author: Paulovich lab**

**Date: 7/26/2016**

## Purpose

The purpose of this document is to describe the liquid chromatography (LC) method for quantitative analysis.

## Scope

This procedure encompasses the setup of the LC and method parameters. It is specific to operation of the Eksigent Ultra system in trap/elute mode.

## Responsibilities

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

## Equipment

- HPLC: UltraNanoLC-2D plus (Eksigent, 950-00061)
- Eksigent AS2 autosampler.

## Materials

- Injection loop: 10 uL peeksil
- Trap column: 300  $\mu\text{m}$  ID x 5 mm LC Packings trap column (Acclaim PepMap 100 C18, 5  $\mu\text{m}$ , 100  $\text{\AA}$ , 160454, Dionex) in a precolumn holder (160431, Dionex)
- Analytical column: 75  $\mu\text{m}$  ID, 360  $\mu\text{m}$  OD IntegraFrit column (IF360-75-50-N-5, New Objective), packed in-house to 10.0 cm with ReproSil-Pur C18-AQ, 3  $\mu\text{m}$  (Dr. Maisch GmbH, Ammerbuch, Germany)
- Water, Optima® LC/MS, suitable for UHPLC-UV (W6-4, Fisher Scientific)

- Acetonitrile (ACN), Optima® LC/MS, suitable for UHPLC-UV (A955-4, Fisher Scientific)
- Formic Acid (FA) (EDM, 11670-1)

## Solutions

Mobile phases. Must be degassed every week:

- Gradient 1, mobile phase A: 0.1% FA in H<sub>2</sub>O
- Gradient 1, mobile phase B: 0.1% FA in 90% ACN

## Procedure

### 1. Autosampler method:

| Step # | Operation       | Value    | Parameter                   | Speed | Height |
|--------|-----------------|----------|-----------------------------|-------|--------|
| 1      | Valve           |          | Injector Load               |       |        |
| 2      | External Events |          | Wait for Grad 1 Ready       |       |        |
| 3      | External Events |          | Wait for Grad 2 Ready       |       |        |
|        | External Events |          | Wait for Loading Pump Ready |       |        |
| 4      | Aspirate        | 11       | Reagent-1                   | 2     | 3      |
| 5      | Wait            | 00:00:02 |                             |       |        |
| 6      | Aspirate        | 0        | Reagent-1                   | 1     | 3      |
| 7      | Aspirate        | 11       | Sample                      | 2     | 2      |
| 8      | Wait            | 00:00:05 |                             |       |        |
| 9      | Aspirate        | 0        | Sample                      | 1     | 2      |
| 10     | Aspirate        | 0        | Reagent-2                   | 1     | 3      |
| 11     | Aspirate        | 2        | Reagent-1                   | 1     | 3      |
| 12     | External Events |          | Start Grad 1                |       |        |
| 13     | External Events |          | Start Grad 2                |       |        |
| 14     | External Events |          | Start Loading Pump          |       |        |
| 15     | Valve           |          | Injector Inject             |       |        |
| 17     | Wait            | 00:05:00 |                             |       |        |
| 18     | Valve           |          | Injector Load               |       |        |
| 19     | Dispense        | 24       | Waste                       | 3     | 0      |
| 20     | Needle Wash     | 100      | Port 1                      |       |        |
| 21     | END             |          |                             |       |        |

### 2. Loading pump method:

- Flow rate ( $\mu\text{L}/\text{min}$ ): See timetable
- Temperature (C): 45
- Run Conditions:
  - Pre-run
  - Flush column for 0.1 minutes using 100% initial flowrate conditions.

## iv. Timetable

| Time (min) | Flow rate (μL/min) | Event |
|------------|--------------------|-------|
| 0          | 20                 |       |
| 10         | 20                 |       |
| 10.1       | 5                  |       |
| 18         | 5                  |       |
| 18.1       | 20                 |       |
| 20         | 20                 |       |

## 3. Gradient 1 method (channel 1):

- i. Flow rate (nL/min): 300
- ii. Temperature (C): 45
- iii. Run Conditions:
  - a. Pre-run
  - b. Flush column for 0.1 minutes using 100% initial flowrate conditions.
- iv. Timetable for column 1 elution:

| Time (min) | % Mobile phase A composition | % Mobile phase B composition | Event        |
|------------|------------------------------|------------------------------|--------------|
| 0          | 98                           | 2                            | Valve Load   |
| 1          | 98                           | 2                            | Valve Inject |
| 2          | 86                           | 14                           |              |
| 12         | 64                           | 36                           |              |
| 13         | 10                           | 90                           |              |
| 14         | 10                           | 90                           |              |
| 15         | 98                           | 2                            |              |
| 20         | 98                           | 2                            |              |