

# STANDARD OPERATING PROCEDURE

**Title: Trypsin Digestion of Cell Lysate, using automated liquid handler**

**Version #: 2 (500 $\mu$ g input)**

**Author: Paulovich lab**

**Date: 8/17/2015**

## Purpose

The purpose of this document is to describe enzymatic digestion of a cell lysate for protein analysis compatible with mass spectrometry.

## Scope

This procedure may be used to reduce, alkylate, and proteolyze samples.

## Responsibilities

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

## Equipment

- EpMotion Automated Pipetting System (Eppendorf).

## Materials

- Trypsin, sequencing grade, (Promega, V5111)
- Urea, (Sigma, U0631)
- 0.5M TCEP, (Pierce, 77720)
- Iodoacetamide (IAM), (Sigma, A3221-10VL)
- 1M Tris (pH8.0), (Sigma, T2694)
- EDTA (Sigma, E7889)
- EGTA (Sigma, E0396)
- Water, HPLC grade (Fisher, W5-1)
- Acetonitrile, HPLC grade (Fisher, A998-1)
- Formic Acid (EDM, 11670-1)

- Sigma Phosphatase Inhibitor Cocktail 1 (Sigma, P2850)
- Sigma Phosphatase Inhibitor Cocktail 2 (Sigma, P5726)
- Deep-well 96-well plate (Eppendorf, 0030 502.248)
- Oasis HLB 96-well plate 30 $\mu$ m (5mg) (Waters, 186000309)
- Square well deep-well V-bottom 96-well plate (Thermo, 95040450).

## Solutions

- Lysis Buffer. **Must be made fresh daily:**
  - 4 Parts 7.5M Urea (see below).
  - 1 Part 5x Lysis Buffer Stock Solution (see below).
  - Add 1% Sigma phosphatase cocktail 1.
  - Add 1% Sigma phosphatase cocktail 2.
  - Mix well.
- 5x Lysis Buffer Stock Solution. May be made in advance and stored at room temp.
  - 12.5mL 1M Tris (pH8.0).
  - 1.0mL 0.5M EDTA.
  - 1.0mL 0.5M EGTA.
  - Add HPLC water to 100mL.
  - Sterilize with 0.22 $\mu$ m filter.
- 7.5 M Urea. **Must be made fresh daily:**
  - Add 4.50g Urea to a 15mL Falcon tube.
  - Add 6mL HPLC water and mix until Urea is in solution.
  - Add HPLC water to a final volume of 10mL.
- 0.2M Tris, pH8.0:
  - 4 parts HPLC water.
  - 1 part 1M Tris, pH8.0.
- 0.5M Iodoacetamide (IAM) Stock Solution. **Prepare immediately before use and keep out of light:**
  - To one 56mg vial of iodoacetamide, add 605 $\mu$ L of 0.2M Tris pH 8.
  - Mix until dissolved.
- Trypsin, Sequence Grade (Promega):
  - 100 $\mu$ g resuspended in 200 $\mu$ L of 0.2M Tris pH 8. (see above).
  - **OR** thaw an aliquot @ 0.525ug/ $\mu$ L.
- Quench:
  - 4 parts HPLC water.
  - 1 part formic acid.
- Heavy peptide stock solutions: Heavy peptides may be obtained from synthetic peptide vendors. Prepare stock solutions of individual heavy peptides at 100-500  $\mu$ M (pmol/ $\mu$ L) in Peptide storage solution (30% acetonitrile / 0.1% formic acid) in sealed screw-top microcentrifuge tubes. Store aliquots at -80  $^{\circ}$ C for 6-12 months.

- Heavy peptide working solution (SIS mix): Generate a master mix of the standard heavy peptides to yield an equal molar mix at 2 pmol/ $\mu$ L in 3% acetonitrile / 0.1% formic acid. Store aliquots at -80 °C for 3 months.

## Procedure

### Preparation of Samples

1. Cell lysate samples were prepared as described in SOP P-02 (Cell line lysis).
2. Adjust the cell lysate to 5mg/mL with lysis buffer.
3. Move 100uL of cell lysate to a deep well plate.
4. Add 6uL of 500mM TCEP.
5. Incubate with mixing at 700rpm for 30 min at 37 °C.
6. Add 14 $\mu$ L of 0.5 M IAM, to yield an IAM concentration of ~40 mM.
7. Alkylate at room temperature for 30 min in the dark.
8. Add 880uL of 200mM Tris, pH 8.0 to each tube to decrease the urea concentration to ~0.6M.
9. Add 20uL trypsin to each digest to achieve a 1:50 enzyme-to-substrate ratio.
10. Incubate with mixing at 700rpm for 2 h at 37 °C.
11. Add 10uL trypsin to each digest to achieve a 1:100 enzyme-to-substrate ratio.
12. Incubate with mixing at 700rpm for 16 h at 37 °C.
13. Add 54 $\mu$ L of 20% FA to each digest to quench the digestion for a final acid concentration of 1%.

### Heavy SIS spike addition

1. The stock heavy SIS mix is at 2uM.
2. Dilute the stock SIS mix to 10 fmol/uL with 3% acetonitrile/0.1% formic acid.
3. 10uL of the diluted spike sample is mixed into each sample.

### Desalting Samples Offline by Positive Pressure

1. Desalt the digested peptides on an Oasis HLB 96 well plate, or equivalent solid phase extraction medium.
  - a. Condition the wells with 4  $\times$  0.4 mL of 0.1% formic acid in 50% acetonitrile.
  - b. Equilibrate the wells using 4  $\times$  0.4 mL of 0.1% formic acid in water.
  - c. Apply the digests to the wells.
  - d. Wash with 4  $\times$  0.4 mL of 0.1% formic acid in water.
  - e. Elute peptides with 3  $\times$  0.4 mL of 0.1% formic acid in 50% acetonitrile. Collect eluates in a 2 mL polypropylene 96 deep well plate.
2. Lyophilize the digests and store frozen (-70 °C) as dried powder until the next step.

## Referenced Documents

SOP P-02 Cell line lysis.pdf