

# STANDARD OPERATING PROCEDURE

**Title: Cell lysate preparation for samples compatible with mass spectrometry**

**Version #: 1**

**Author: Paulovich Lab**

**Date: 7/26/16**

## Purpose

The purpose of this document is to describe preparation of a cell lysate for protein analysis compatible with mass spectrometry.

## Scope

This procedure may be used to make a cell lysate from cell lines.

## Responsibilities

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

## Equipment

- Refrigerated micro-centrifuge
- Centrifuge for use with tubes and 96 well plates
- Coulter Counter or hemocytometer
- Sonic dismembrator, (model 550, ThermoFisher Scientific)
- Lyophilizer

## Materials

- UREA: SIGMA ULTRA CAT# U0631
- 1 M TRIS (PH8.0): SIGMA CAT# T2194
- EDTA: SIGMA CAT# E7889
- EGTA: SIGMA CAT# E0396
- HPLC WATER: FISHER CAT# W6-4

- SIGMA PHOSPHATASE COCKTAILS (P2850, P5726)

## Reagents and Solutions

1. 1×DPBS (Gibco, ThermoFisher Scientific) (138 mM NaCl, 2.7 mM KCl, 8 mM Na<sub>2</sub>HPO<sub>4</sub>·2H<sub>2</sub>O, 1.5 mM KH<sub>2</sub>PO<sub>4</sub>).
2. 7.5 M Urea: Add 6 mL HPLC-grade water to a 15 mL Falcon conical tube. Add 4.50 g urea and mix until the urea is in solution. Bring final volume to 10 mL with HPLC-grade water. Must be made fresh daily. Store at room temperature.
3. 5×Lysis buffer stock solution: To Make 100 mL: Add 12.5 mL 1 M Tris (pH 8.0), 1.0 mL 0.5 M EDTA, and 1.0 mL 0.5 M EGTA. Bring to 100 mL with HPLC-grade water. Filter with 0.22 μm filter. Store at room temperature for 6 months.
4. Lysis buffer: Add 4 Parts 7.5 M Urea, 1 Part 5×Lysis buffer stock solution, 1% Sigma phosphatase cocktail 2, 1% Sigma phosphatase cocktail 3, and 1% Sigma Protease Inhibitor to a final volume of 100 mL. Mix well. Must be made fresh daily, store on ice.
5. 0.2 M Tris, pH 8.0. Add 4 parts HPLC-grade water and 1 part 1 M Tris, pH 8.0.

## Procedure

1. Preparation:
  - a. Turn on bench-top centrifuge and cool to 4°C (check that appropriate adaptors are in the buckets).
  - b. Turn on refrigerated micro-centrifuge and cool to 4°C.
  - c. Turn on Coulter Counter and prime aperture.
  - d. Thaw Phosphatase inhibitors.
  - e. Label and pre-cool 50 mL tubes.
  - f. Label and pre-cool micro-centrifuge tubes.
  - g. Label and pre-cool Cryo-vials.
  - h. Make fresh urea lysis buffer- see solution section above.
2. Remove media and rinse cells with 10 mL DPBS.
3. Harvest ~ two T-175 flasks by Trypsinization (6 mL / T175 flask or 2 mL / plate)
4. Incubate cells at room temperature with occasional mixing until the cells lift from culture surface as seen under the microscope.
5. Add 20 mL / flask (or 8 mL / plate) Growth Media (10% FBS) to swamp out the Trypsin
6. If needed, use a cell lifter to remove all cells from culture surface.
7. Transfer cells to pre-cooled 50 mL tubes.
8. Spin cells 180 x g / 8 min. / 4°C, dump supernatant.
9. Resuspend and pool cells from the same cell line in 10 mL ice-cold DPBS, remove 50 uL for cell counting by Hemocytometer.
10. Add ice-cold DPBS to 50 mL
11. Spin cells 180 x g / 8 min. / 4°C.
12. Count cell aliquot while cells are spinning. Calculate total cell count.

13. Dump supernatant and resuspend cells in 50 mL of ice-cold DPBS.
14. Spin cells 180 x g / 8 min. / 4°C.
15. Dump supernatant and invert tubes on paper towel ~ 20 sec., tap tubes on paper towels and remove remaining supernatant w/ Kim-Wipe.
16. Prepare Lysate buffer (see Solution section above)
  - a. Lyse at  $0.5 \times 10^8$  cells / mL
  - b. Lysate buffer needs to be made fresh just before use.
  - c. Add phosphatase inhibitors to lysis buffer at 1% just before use.
17. Place tube with cell pellet on ice and add lysis buffer to a final concentration of  $0.5 \times 10^8$  cells / mL.
18. Gently resuspend cells in lysis buffer by dragging tube along a microfuge tube rack- do not pipette.
19. Sonicate cells 2 x 10 sec. (550 Sonic Dimembrator, Fisher Scientific; knob set to (5)
  - a. Wipe down probe with water and ethanol between samples.
  - b. Place lysate on ice for ~10 sec. between sonications.
20. Transfer lysate by pipette tip to micro-centrifuge tube, vortex 15 sec., ice 10 min., vortex 15 sec.
21. Micro-centrifuge / 20k x g (14K RPM or full speed) / 10 min. / 4°C.
22. Transfer supernatant to 1.0 mL cryo-vial (Nunc Cat# 377267).
  - a. Note: if storing aliquots of a lysate, first transfer the lysate to a fresh micro-centrifuge tube to ensure homogeneous mixing of the lysate before aliquoting.
23. Determine initial protein concentration by Bradford or BCA assay (dilute lysate 1:25 in H<sub>2</sub>O).
24. Store lysates in liq. N<sub>2</sub>.
25. When all lysates are available for pooling (if necessary), thaw lysates on ice.
26. Determine protein concentration by BCA.
27. Pool lysates as desired.

## Referenced Documents

n/a