

STANDARD OPERATING PROCEDURE

Title: Peptide immunoaffinity enrichment for the characterization of immuno-MRM assays using response curves

Version #: 1

Author: Paulovich Lab

Date: 9/18/2014

Purpose

The purpose of this document is to describe the immunoaffinity enrichment of peptides as part of the characterization of immuno-MRM assays using response curves.

Scope

This procedure may be used to generate immunoaffinity-enriched peptide response curve samples to be analyzed by MRM. It includes a description of i) generating the response curve samples, ii) setting up the immunoaffinity captures using the samples and a mix of antibodies that had previously been cross-linked to magnetic beads, and iii) the automated washing and elution of the target peptides using a robotic magnetic bead handling system. Last, the run order of the samples on the LC-MRM-MS system is given.

Responsibilities

It is the responsibility of person(s) performing this procedure to be familiar with laboratory safety procedures. The interpretation of results must be done by a person trained in the procedure and familiar with such interpretation.

Equipment

- KingFisher 96 magnetic bead handling platform (5400500, Thermo Scientific), equipped with a PCR magnet head
- LabQuake tube rotator (415110, Barnstead)
- Dynabeads® MPC®-S (Magnetic Particle Concentrator, A13346, Life Technologies) or DynaMag™-15 Magnet (12301D, Life Technologies)
- DynaMag™-96 Side Skirted Magnet (12027, Life Technologies)
- Centrifuge for use with 96 well plates (5810R, Eppendorf)

Materials

- 1.7 mL microcentrifuge tubes (24-282, Genesee Scientific)
- 15 mL conical centrifuge tubes (352196, Corning Life Sciences)
- MColorpHast® pH Test Strips, 5-10 (1.09533.0001, EMD Millipore)
- KingFisher 96 KF microplate (200 μ L) (97002540, Thermo Scientific)
- KingFisher 96 tip comb for PCR magnets (97002514, Thermo Scientific)
- Hard-Shell® Low-Profile Thin-Wall 96-Well Skirted PCR Plates (HSP-9601, BioRad)
- Chemically resistant sealing foil (T-3025-8B, BioExpress)
- ThermalSeal® PCR Sealing Films (12-168, Genesee Scientific)
- 1 M Trizma hydrochloride buffer, pH 8.0 (Tris, T2694, Sigma)
- 10X Phosphate buffered saline (PBS) (BP39920, Fisher Scientific)
- CHAPS detergent (28300, Thermo Scientific)
- Water, Optima® LC/MS, suitable for UHPLC-UV (W6-4, Fisher Scientific)
- Acetic acid, ACS reagent, \geq 99.7% (242853, Sigma)
- Acetonitrile, Optima® LC/MS, suitable for UHPLC-UV (A955-4, Fisher Scientific)
- Sodium azide (NaN_3) (S8032, Sigma)
- Lima bean trypsin inhibitor (LS002829, Worthington Biochemical Corporation)
- Bovine 6 protein tryptic digest equal molar mix (PTD/00006/63, Bruker)

Reagents

- Light and stable isotope-labeled heavy synthetic peptides for MRM analyses were purchased at $>$ 95% purity as measured by HPLC. Cysteine residues of the tryptic light and heavy peptides were carbamidomethylated, and the heavy peptides' C-terminal lysine (K) or arginine (R) residues were uniformly ^{13}C and ^{15}N labeled, yielding a mass increase of 8 Da (for K) or 10 Da (for R). For peptides at the C-terminus of proteins, other amino acids (such as valine) can be isotopically labeled. Amino acid analysis (AAA) was performed to determine peptide concentrations. Individual peptide stock aliquots were stored in 30% acetonitrile-0.1% formic acid at 100 pmol/ μ L at -80°C .
- Rabbit-anti-human monoclonal antibodies (mAbs) were custom-generated to the tryptic and proteotypic peptide targets and covalently coupled to magnetic Protein G beads. For the coupling procedure, see the SOP titled "Antibody coupling and immobilization to Protein G magnetic beads using a 96-well plate KingFisher automated magnetic bead handling platform."
- Trypsinized human plasma was prepared according to the SOP titled "Trypsin Digestion of Human Plasma."

Solutions

- Prepare approximately 30 μ L of a 0.1 mg/ μ L solution of lima bean trypsin inhibitor in Optima® water.

- Prepare a 5% (w/v) CHAPS solution in Optima® water by dissolving 1.0 g of CHAPS powder in 20 mL of Optima® water. It will take some time to dissolve the powder completely. (CHAPS aids in the bead transfer from plate to plate on the KingFisher.)
- Prepare 100 mL of a 1X PBS-0.01% CHAPS solution, pH 7 by adding 10.0 mL of 10X PBS and 200 μ L of 5% CHAPS to Optima® water, for a total of 100 mL.
- Prepare 100 mL of a 1/10X PBS-0.01% CHAPS solution, pH 7 by adding 1.0 mL of 10X PBS and 200 μ L of 5% CHAPS to Optima® water, for a total of 100 mL.
- Prepare 100 mL of 3% acetonitrile-5% acetic acid by adding 3.0 mL of Optima® acetonitrile and 5.0 mL of glacial acetic acid to Optima® water, for a total of 100 mL.
- Prepare a 10% (w/v) NaN_3 solution in water by dissolving e.g. 1.0 g of NaN_3 powder in 1.0 mL of Optima® water.
- Prepare 50 mL of a 1X PBS-0.03% CHAPS-0.1% NaN_3 solution by adding 5.0 mL of 10X PBS, 300 μ L of 5% CHAPS, and 500 μ L of 10% NaN_3 to Optima® water, for a total of 50 mL.
- Prepare a light peptide master mix with a 2 pmol/ μ L concentration per peptide from the individual light peptide stocks (at 100 pmol/ μ L, at -80 °C). For 40 individual peptides, for example, this would be done by adding 30 μ L of Optima® water to a 1.7 mL microcentrifuge tube, and then adding 3 μ L of each peptide stock to the tube. Aliquot this master mix into 1.7 mL microcentrifuge tubes, with 10 μ L of the master mix in each aliquot, and store at -80 °C until use.
- Prepare, aliquot, and store a heavy peptide master mix with a 2 pmol/ μ L concentration per peptide as was done for the light peptide master mix.

Procedure

The following describes the experiments for a 9-point response curve that has varying concentrations of the heavy (analyte) peptides and a constant concentration of the light (internal standard) peptides in a human plasma digest matrix, including a 0 fmol (blank) sample that does not contain heavy peptides. The 0 fmol sample is prepared for 5 replicate immuno-MRM experiments, and the other samples are prepared for triplicate experiments per curve point. The light peptides were not chosen as analytes since the endogenous detection of some of the peptides in human plasma would not have made it possible to determine the limits of detection and quantification.

The experiments encompass

- the preparation of a monoclonal antibody master mix (with each antibody having previously been covalently coupled (cross-linked) to magnetic Protein G beads, see the SOP titled “Antibody coupling and immobilization to Protein G magnetic beads using a 96-well plate KingFisher automated magnetic bead handling platform”),

- the preparation of the different curve point samples (with light and heavy peptides spiked into a diluted human plasma digest matrix, with 10 μL of original plasma used per immunoaffinity enrichment (capture) experiment),
- the immunoaffinity enrichment of the peptides from the plasma digest samples, and
- the automated washing of the beads+antibodies+peptides and the elution of the peptides from the antibodies using a KingFisher magnetic bead handling platform.

Each capture sample contains 200 μL of a mix of peptides and plasma digest, plus 50 μL of the antibody master mix. These volumes have been convenient to use for this procedure, yet could be changed to e.g. 150 μL and 25 μL , respectively. Of importance is, however, that the *concentrations* used for plotting the response curves in the end are defined as X fmol of heavy peptide per 10 μL of plasma, or X/10 fmol of peptide per 1 μL of plasma.

Preparation of the monoclonal antibody (mAb) master mix

1. Calculate the *total volume* of the antibody master mix needed based on the 50 μL of antibody-beads master mix wanted per capture and the number of captures for the response curves (allowing for a couple of extra captures for good measure).
2. Calculate the volume of *each individual* cross-linked antibody-to-beads stock needed for all captures, given an antibody concentration of 0.5 μg antibody/ μL of the individual antibody-beads suspension (see the SOP titled “Antibody coupling and immobilization to Protein G magnetic beads using a 96-well plate KingFisher automated magnetic bead handling platform”) and the number of captures (plus 2 extra captures for good measure).

In the end, depending on which volume from these two calculations is larger, excess supernatant (1X PBS-0.03% CHAPS-0.1% NaN_3) will either need to be removed from the master mix, or extra 1X PBS-0.03% CHAPS-0.1% NaN_3 can be added to the master mix to arrive at the 50 μL of antibody master mix volume for each capture. (Alternatively, since the 50 μL per capture volume is rather arbitrary, 1X PBS-0.03% CHAPS-0.1% NaN_3 does not need to be added and a lower volume of the antibody master mix can be used per capture.)

3. Add the calculated volume of each antibody-beads suspension stock to either a 1.7 mL microcentrifuge tube or a 15 mL conical tube.
4. If excess supernatant needs to be removed, put the tube onto a magnetic particle concentrator for 2 minutes and remove the excess supernatant by pipetting.
5. Vortex the tube and spin it down briefly, then keep it on ice until use.

Preparation of the response curve samples

The levels of heavy peptides that were tested in these particular response curves were 5000, 1000, 200, 40, 8, 2, 0.5, 0.125, and 0 fmol per 10 μL of plasma digest. These levels were generated using serial dilutions, with the dilutions from 5000 to 8 fmol/10 μL plasma being covered by 1:5 serial dilutions, and from 8 to 0.125 fmol/10 μL plasma by 1:4 serial

dilutions (smaller intervals were chosen for the lower levels since this would allow more robust characterizations of the limits of quantification). The synthetic light peptides were spiked at a constant level of 100 fmol per 10 μL of plasma digest. Considering a) the 200 μL volume of peptides-in-plasma wanted per capture sample, b) the number of replicates per curve point, c) wanting to be able to keep the same pipettor volume setting during the serial dilution, and d) wanting comfortably enough volume at each step for pipetting purposes, the following procedure was followed to prepare the response curve samples.

6. Thaw a tube containing 1 mL of human plasma digest (see SOP titled "Trypsin Digestion of Human Plasma") on ice and add 15 μL of lima bean trypsin inhibitor solution at 0.1 mg/ μL to yield a 1:2 trypsin inhibitor:trypsin ratio (w/w) (to avoid potential digestion of the antibodies during the enrichment).
7. Prepare a light peptide-in-plasma digest mix as follows (this yields a total volume of 8600 μL):
 - 7.1. Add 7738 μL of 1X PBS-0.01% CHAPS to a 15 mL conical centrifuge tube;
 - 7.2. add 430 μL of human plasma digest (equivalent to 430 μL of original, undigested plasma);
 - 7.3. add 430 μL of 1 M Tris, pH 8.0 (to ensure a pH between 7.0 and 8.0 since the plasma digest can be acidic after digestion and desalting; the antibodies that are added later might degrade under acidic conditions);
 - 7.4. add 2.15 μL of the 2 pmol/ μL synthetic light peptide master mix (yielding a synthetic light peptide concentration of 100 fmol/10 μL of plasma).
 - 7.5. vortex and spin down the sample briefly.
8. Aliquot the light peptide-in-plasma digest mix, add heavy peptides, and perform the serial dilutions as follows:
 - 8.1. Aliquot the light peptide-in-plasma mix into eight 1.7 mL microcentrifuge tubes:
 - a) 1100 μL into 1 tube (label the tube with 5000 (fmol/10 μL plasma))
 - b) 880 μL into each of 4 tubes (label the tubes with 1000, 200, 40, and 8 (fmol/10 μL plasma))
 - c) 660 μL into each of 3 tubes (label the tubes with 2, 0.5, and 0.125 (fmol/10 μL plasma)).
 - 8.2. Place the tube with the remaining light peptide-in-plasma digest mix on ice (5 x 200 μL of this mix will later be transferred to 5 wells of a 96 well plate, representing the 0 fmol/10 μL plasma replicate samples).
 - 8.3. Take 13.75 μL out of the first (5000) tube and discard, and add 13.75 μL of the 2 pmol/ μL heavy peptide master mix to that tube; mix by centrifuging and spin down briefly. Since the heavy peptide mix contained formic acid, measure the pH of the peptide+plasma sample using e.g. 2 μL spotted onto a pH 5.0-10.0 pH test strip to ensure a pH above 7.0. The heavy peptide concentration in this sample is 5000 fmol/10 μL of plasma.

8.4. Perform the 1:5 and 1:4 serial dilutions using 220 μL from the first tube and diluting first in the tubes containing 880 μL of the light peptide-in-plasma mixes, and then into the tubes containing 660 μL . Ensure good mixing at each dilution step.

Overnight immunoaffinity enrichment (capture) of the peptides from the plasma digest samples

9. Add 3x 200 μL of each concentration's response curve sample to 3 wells of a KingFisher 96 KF microplate (cat. # 97002540; the product description calls it a 200 μL -per-well plate, but more than 200 μL of sample can be added if a KingFisher 96 tip comb for PCR magnets is used; here, the total volume will be 250 μL once the beads+antibodies are added).
10. Add 5x 200 μL of the light peptide-in-plasma digest mix to 5 wells of the microplate (representing the 0 fmol/10 μL plasma replicate samples).
11. Vortex the tube containing the antibody master mix until all beads are well suspended. Add 50 μL of the master mix to each well using a single-channel pipette, re-capping the tube after every $\sim 5^{\text{th}}$ addition and briefly vortexing it to continue to ensure a well-mixed bead suspension. Alternatively, the beads can also be added using a multi-channel pipette and a pipette tray.
12. If beads have settled to the bottom of the plate wells at the end of all additions, they can be briefly suspended again using a multi-channel pipette.
13. Immediately cover the plate using chemically resistant sealing foil, making sure that each well is completely sealed so no liquid can escape from any of the wells.
14. Quickly secure the plate onto a LabQuake rotator using lab tape, and move the rotator to a fridge or cold room and rotate the plate overnight at 4 $^{\circ}\text{C}$.

Automated washing of the beads+antibodies(Abs)+peptides and elution of the peptides from the antibodies using a KingFisher magnetic bead handling platform

The following is an overview of the types of plates and contents that will be used with the KingFisher platform. Plate 1 should be left rotating on the LabQuake until the other plates have been prepared.

Plate #	Plate type	KingFisher step	Contents	Volume per well (μL)
1	KingFisher 96 KF microplate	Bead+Ab+peptide mixing and pick-up	Response curve samples	250
2	KingFisher 96 KF microplate	Wash 1	1X PBS-0.01% CHAPS	250
3	KingFisher 96 KF microplate	Wash 2	1X PBS-0.01% CHAPS	250
4	KingFisher 96 KF microplate	Wash 3	1/10X PBS-0.01% CHAPS	250
5	Hard-Shell® Low-Profile Thin-Wall 96-Well Skirted PCR Plates	Elution	3% acetonitrile-5% acetic acid	26
6	KingFisher 96 KF microplate	Bead+Ab collection	1X PBS-0.03% CHAPS-0.1% sodium azide	100
7	KingFisher 96 KF microplate	Tip comb	N/A	N/A

The following shows the details of the KingFisher method employed for washing and eluting the samples. No precollection was employed at the beginning of any step, the samples were not heated during mixing, and no postmixing was employed at the end of any step. The method takes about 25 minutes.

Plate #	KingFisher step	Beginning of step release time (hr:min:sec), speed	Mixing time (hr:min:sec), speed	End of step collect count	End of step collect time (sec)
7	Tip comb; Pick-up	N/A	N/A	N/A	N/A
1	Bead+Ab+peptide mixing and pick-up	00:00:10, Slow	00:05:00, Bottom slow	5	1.5
2	Wash 1	00:00:10, Slow	00:01:30, Slow	5	1.5
3	Wash 2	00:00:10, Slow	00:01:30, Slow	5	1.5
4	Wash 3	00:00:10, Slow	00:01:30, Slow	5	1.5
5	Elution	00:00:10, Bottom slow	00:05:00, Bottom slow	10	1.5
6	Bead+Ab collection	00:00:20, Fast	N/A	N/A	N/A
7	Tip comb; Leave	N/A	N/A	N/A	N/A

15. Prepare Plate 5, cover it with a ThermalSeal® PCR Sealing Film, spin it down briefly, and place it on ice.
16. Prepare Plates 2, 3, 4, 6, and 7 and place them into the KingFisher.
17. Take the plate that contains the response curve samples (Plate 1) off the LabQuake and spin it down at 1000 rpm for 15 sec. This will spin down any liquid that might have contacted the foil cover.
18. Carefully take the foil cover off Plate 1, then resuspend the beads that have been compacted at the bottom of the wells using a multichannel pipette.
19. Right after the resuspension, place Plate 1 into the KingFisher and start the KingFisher method (this is done without delay to avoid having the beads settle again, which can cause them to not be transferred completely to Plate 2 (a small amount of beads might still be left behind in Plate 1, but this should be negligible compared to the total amount of beads)).

20. Once the KingFisher starts the first wash step, take the PCR Sealing Film off Plate 5 and place the plate into the KingFisher.
21. Once the method has finished, take Plate 5 out of the KingFisher, cover the plate with PCR Sealing Film, and centrifuge it briefly.
22. Place Plate 5 onto a DynaMag™-96 Side Skirted Magnet for 5 minutes, then split the eluate samples into 2 new PCR plates, with 13 μ L of each sample per plate (11 μ L will be used per LC-MRM-MS run for a full-loop injection using a 10 μ L sample loop). Spin down both plates and freeze them at -80 °C until use. One plate is a back-up plate in case a second injection should be needed; the back-up plate should not be thawed unless necessary.
23. The beads+Abs that were collected in plate 6 may be saved and stored at 4 °C.

Sample run order on the LC-MRM-MS system

The descriptions of LC and MRM system set-ups are given in the SOPs titled “Liquid Chromatography using an Eksigent NanoLC-2D, an Eksigent AS1 autosampler, two Agilent 10-port valves, a C18 trap column, and a ReproSil C18 analytical column for immuno-MRM experiments” and “MRM mass spectrometry for the analysis of immuno-MRM assay samples using a 5500 QTRAP,” respectively.

The response curve samples are run in increasing order of heavy peptide level, with one injection per capture replicate. Before that, a system suitability standard peptide mix (unrelated peptides, e.g. a bovine 6 protein tryptic digest equal molar mix) is run to ascertain that the system is behaving as expected (see also Abbatiello, S. E., *et al.*¹). An injection of the synthetic light peptide standards in buffer is then performed to determine the approximate retention times of the peptides, with which a scheduled LC-MRM method is generated using the Skyline software². This is followed by the samples, with a buffer blank run (injecting 10 μ L of just 3% acetonitrile-5% acetic acid) in between each sample run. The table below illustrates the run order. (Due to some peptides showing considerable carry-over from run to run, it is not practical to randomize the run order with respect to heavy peptide spike level.)

Sample	Notes
System suitability standard peptide mix, run 1	E.g. a 6 bovine protein tryptic digest equal molar mix
System suitability standard peptide mix, run 2	Checking for consistency of the peptide profile; a third run of this mix might be needed if the system had been idle and the first run was used to equilibrate the system.
Synthetic light peptide mix in buffer	For retention time scheduling purposes.
Wash	10 μ L of 50% methanol-40% acetonitrile-10% Optima® water.
Buffer blank	10 μ L of 3% acetonitrile-5% acetic acid
0 fmol / 10 μ L plasma sample, repeat 1	Adjust ret. time scheduling based on this run if needed.
Buffer blank	
Buffer blank	This run is used to characterize the light peptide carry-over signal that remains after the first buffer blank run; ideally, the light peptide intensity or peak area remaining in this run should be no more than 2% of the signal in the 0 fmol / 10 μ L plasma sample, repeat 1.
0 fmol / 10 μ L plasma sample, repeat 2	
Buffer blank	
0 fmol / 10 μ L plasma sample, repeat 3	
Buffer blank	
0.125 fmol / 10 μ L plasma sample, repeat 1	
Buffer blank	
0.5 fmol / 10 μ L plasma sample, repeat 1	
Buffer blank	
2 fmol / 10 μ L plasma sample, repeat 1	
Buffer blank	
etc.	
5000 fmol / 10 μ L plasma sample, repeat 1	
Buffer blanks and washes	As many runs as needed until the heavy peptide signal has returned to baseline (i.e., no heavy peptide is detected).
0 fmol / 10 μ L plasma sample, repeat 4	
Buffer blank	
0.125 fmol / 10 μ L plasma sample, repeat 2	
Buffer blank	
0.5 fmol / 10 μ L plasma sample, repeat 2	
etc.	

Referenced Documents

SOP “Antibody coupling and immobilization to Protein G magnetic beads using a 96-well plate KingFisher automated magnetic bead handling platform”

SOP “Trypsin Digestion of Human Plasma”

SOP “Peptide immunoaffinity enrichment for the characterization of immuno-MRM assays using response curves”

SOP “Liquid Chromatography using an Eksigent NanoLC-2D, an Eksigent AS1 autosampler, two Agilent 10-port valves, a C18 trap column, and a ReproSil C18 analytical column for immuno-MRM experiments”

1. Abbatiello, S. E., Mani, D. R., Schilling, B., Maclean, B., Zimmerman, L. J., Feng, X., Cusack, M. P., Sedransk, N., Hall, S. C., Addona, T., Allen, S., Dodder, N. G., Ghosh, M., Held, J. M., Hedrick, V., Inerowicz, H. D., Jackson, A., Keshishian, H., Kim, J. W., Lyssand, J. S., Riley, C. P., Rudnick, P., Sadowski, P., Shaddock, K., Smith, D., Tomazela, D., Wahlander, A., Waldemarson, S., Whitwell, C. A., You, J., Zhang, S., Kinsinger, C. R., Mesri, M., Rodriguez, H., Borchers, C. H., Buck, C., Fisher, S. J., Gibson, B. W., Liebler, D., MacCoss, M., Neubert, T. A., Paulovich, A., Regnier, F., Skates, S. J., Tempst, P., Wang, M., and Carr, S. A. (2013) Design, implementation and multisite evaluation of a system suitability protocol for the quantitative assessment of instrument performance in liquid chromatography-multiple reaction monitoring-MS (LC-MRM-MS). *Mol Cell Proteomics* 12, 2623-2639.
2. MacLean, B., Tomazela, D. M., Shulman, N., Chambers, M., Finney, G. L., Frewen, B., Kern, R., Tabb, D. L., Liebler, D. C., and MacCoss, M. J. (2010) Skyline: an open source document editor for creating and analyzing targeted proteomics experiments. *Bioinformatics* 26, 966-968.